


# ROVI REQUEST FOR QUOTE

Please complete entire form



## PROVIDER INFORMATION

<b>VA</b>   	U.S. Department of Veterans Affairs	Billing Address: RECIPIENT NAME:	Shipping Address: RECIPIENT NAME:	Same as billing
Ordering Clinic Name:	STREET:	STREET:		
Ordering Clinician:	CITY:	STATE:	CITY:	STATE:
Client Reference / Order Reference Number:	ZIP CODE:	COUNTRY:	ZIP CODE:	COUNTRY:

## PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or a retrofit to an existing wheelchair?  New Chair Order  Retrofit Existing Chair

**Select LUCI Model:** VA PRICE: \$8,445.00 USD

*Note: We do not support swing-away footrests, the extra-large flip-up foot platform, or high speed motors. If you intend to use specialty controls please contact LUCI to confirm compatibility.*

ROVI X3 - LUCI Model 111

**Select Seat-to-Floor Height:**

*Note: Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. LUCI is not compatible with 17 inch seat to floor height as some seating may interfere with LUCI and/or limit seat travel.*

18 inches      19 inches      20 inches

Any custom chair modifications? If so, explain.

**Select Drive Method:**

*This will ensure you receive the correct brackets*

- R-Net LED Remote Joystick
- R-Net LED Remote Joystick for Lights
- R-Net Color Remote Joystick with Monojacks
- R-Net Color Remote Joystick with Monojacks for Lights
- R-Net Super Joystick
- Omni2: ALT Drive Method

*Note: We do not support Latch drive or VR2 devices*

If you are using a configuration not found on this form, please call **615-813-5824, ext 1** before ordering.

Send completed form to [orders@luci.com](mailto:orders@luci.com)