


ROVI REQUEST FOR QUOTE

Please complete entire form



PROVIDER INFORMATION

VA |  U.S. Department of Veterans Affairs

Billing Address: _____ Shipping Address: *Same as billing*

RECIPIENT NAME: _____ RECIPIENT NAME: _____

Ordering Clinic Name: _____ STREET: _____ STREET: _____

Ordering Clinician: _____ CITY: _____ STATE: _____ CITY: _____ STATE: _____

Client Reference / Order Reference Number: _____ ZIP CODE: _____ COUNTRY: _____ ZIP CODE: _____ COUNTRY: _____

PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or to retrofit an existing wheelchair? New Chair Order Retrofit Existing Chair

Note: For New Chairs, on the ROVI order form select the LUCI Interface Ready option. This will ensure the chair arrives with pre-installed brackets for easy installation. Brackets can also be installed when retrofitting LUCI to an existing chair.

Select LUCI Model: VA PRICE: \$8,445.00 USD

Note: We do not support swing-away footrests, the extra-large flip-up foot platform, or high speed motors. If you intend to use specialty controls please contact LUCI to confirm compatibility.

ROVI X3 - LUCI Model 111

ROVI A3 - LUCI Model 112

Select Seat-to-Floor Height:

Note: Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. LUCI is not compatible with 17 inch seat to floor height as some seating may interfere with LUCI and/or limit seat travel.

18 inches

19 inches

20 inches

Any custom chair modifications? If so, explain.

Select Drive Method:

This will ensure you receive the correct brackets

R-Net LED Remote Joystick

R-Net LED Remote Joystick for Lights

R-Net Color Remote Joystick with Monojacks

R-Net Color Remote Joystick with Monojacks for Lights

R-Net Advanced Joystick

Omni2: ALT Drive Method

Note: We do not support Latch drive or VR2 devices

If you are using a configuration not found on this form, please call **615-813-5824, ext 1** before ordering.

Send completed form to orders@luci.com

ROVI-0523-V