

# ROVI REQUEST FOR QUOTE

Please complete entire form



## PROVIDER INFORMATION

Provider Company Name:	Billing Address:	Shipping Address:	<i>Same as billing</i>
	RECIPIENT NAME:	RECIPIENT NAME:	
Provider Contact Name:	STREET:	STREET:	
Referring Therapist, Clinic, and Contact Info:	CITY:	STATE:	CITY: STATE:
Client Reference / Order Reference Number:	ZIP CODE:	COUNTRY:	ZIP CODE: COUNTRY:

## PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or to retrofit an existing wheelchair?  New Chair Order  Retrofit Existing Chair  
Note: For New Chairs, on the ROVI order form select the LUCI Interface Ready option. This will ensure the chair arrives with pre-installed brackets for easy installation. Brackets can also be installed in the field when retrofitting LUCI to an existing chair.

Select LUCI Model: [MSRP: \\$9854.00 USD](#)  
Note: We do not support swing-away footrests, the extra-large flip-up foot platform, or high speed motors. If you intend to use specialty controls please contact LUCI to confirm compatibility.

- ROVI X3 - LUCI Model 111
- ROVI A3 - LUCI Model 112

Select Seat-to-Floor Height:  
Note: Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. LUCI is not compatible with 17 inch seat to floor height as some seating may interfere with LUCI and/or limit seat travel.

- 18 inches
- 19 inches
- 20 inches

Insurance (K0108 – Miscellaneous) or Private Pay?

- Insurance Carrier
- Private Pay

## REQUEST DELIVERY SUPPORT

LUCI is a finely-tuned, highly customizable system and works best when it is properly setup for the power wheelchair user.

Let LUCI be your guide during delivery and setup and schedule a 30-minute remote support session with us!

Select Drive Method:  
This will ensure you receive the correct brackets

- R-Net LED Remote Joystick
- R-Net LED Remote Joystick for Lights
- R-Net Color Remote Joystick with Monojacks
- R-Net Color Remote Joystick with Monojacks for Lights
- R-Net Advanced Joystick
- Omni2: ALT Drive Method

Note: We do not support Latch drive or VR2 devices

Any custom chair modifications? If so, explain.

Scan here to schedule your session:  
Delivery Contact Name:

Wheelchair User Name:



<https://calendly.com/d/2v7-77p-s64/luci-delivery-support>

Send completed form to [orders@luci.com](mailto:orders@luci.com)

If you are using a configuration not found on this form, please call 615-813-5824, ext 1 before ordering.

ROVI-1023-C