

ROVI REQUEST FOR QUOTE

Please complete entire form



PROVIDER INFORMATION

Provider Company Name:	Billing Address:	Shipping Address:	<i>Same as billing</i>
	RECIPIENT NAME:	RECIPIENT NAME:	
Provider Contact Name:	STREET:	STREET:	
Referring Therapist, Clinic, and Contact Info:	CITY:	STATE:	CITY: STATE:
Client Reference / Order Reference Number:	ZIP CODE:	COUNTRY:	ZIP CODE: COUNTRY:

PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or to retrofit an existing wheelchair? New Chair Order Retrofit Existing Chair

Note: For New Chairs, on the ROVI accessory order form select the LUCI Interface Ready option. This will ensure the chair arrives with pre-installed brackets for easy installation. Brackets can also be installed in the field when retrofitting LUCI to an existing chair.

Select LUCI Model: [MSRP: \\$9854.00 USD](#)

Note: We do not support swing-away footrests, the extra-large flip-up foot platform, lateral tilt, or high speed motors. If you intend to use specialty controls please contact LUCI to confirm compatibility.

ROVI X3 - LUCI Model 111

ROVI A3 - LUCI Model 112

ROVI A3 Stander - LUCI Model 112

Select Seat-to-Floor Height:

Note: Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. LUCI is not compatible with 17 inch seat to floor height as some seating may interfere with LUCI and/or limit seat travel.

18 inches 19 inches 20 inches

Insurance (K0108 – Miscellaneous) or Private Pay?

Insurance Carrier

Private Pay

REQUEST DELIVERY SUPPORT

LUCI is a finely-tuned, highly customizable system and works best when it is properly setup for the power wheelchair user.

Let LUCI be your guide during delivery and setup and schedule a 30-minute remote support session with us!

Select Drive Method:

This will ensure you receive the correct brackets

R-Net LED Remote Joystick

R-Net LED Remote Joystick for Lights

R-Net Color Remote Joystick with Monojacks

R-Net Color Remote Joystick with Monojacks for Lights

R-Net Advanced Joystick

Omni2: ALT Drive Method

Note: We do not support Latch drive or VR2 devices

Any custom chair modifications? If so, explain.

Scan here to schedule your session:

Delivery Contact Name:

Wheelchair User Name:



<https://calendly.com/d/2v7-77p-s64/luci-delivery-support>

Send completed form to orders@luci.com

If you are using a configuration not found on this form, please call 615-813-5824, ext 1 before ordering.

ROVI-0124-C