



# WHEELCHAIR SMARTER.

AWARD-WINNING SMART WHEELCHAIR TECHNOLOGY  
RE-IMAGINING MOBILITY AT EVERY TURN.

**LUCI is the first and only smart technology accessory for power wheelchairs, designed to give riders a safer and more inclusive experience.**


LUCI prevents power wheelchair tips, collisions and falls before they happen, thanks to a platform which helps the wheelchair “see” its environment and respond to steering input from the user. The result is safe, confident independence.




Collision Avoidance



Drop-off and Tip Protection



Cloud-based Communication and Alerts



Enhanced Mobility Platform



Home Assistant Compatibility



AWARDED BY

TIME

POPULAR SCIENCE





## LUCI'S 360° VIEW

**LUCI adds hardware and software to riders' new and existing chairs, offering unparalleled stability, security and connectivity.**

Our patented system combines data from sensor fusion, infrared, ultrasonics and radar into a single view of the world, enabling never-before-seen possibilities for power wheelchair riders.

LUCI is currently available for some of the world's top-selling power wheelchairs, including:

- Permobil M3
- Permobil M5
- Permobil F3
- Permobil F5
- Quickie Q700 M
- Quickie Q500 M
- Quickie Q300 M
- ROVI X3
- ROVI A3 Stander
- Amy Systems R3

SMART WHEELCHAIR TECHNOLOGY FOR  
THE MOST FEARLESS PEOPLE ON EARTH.

## FUELED BY COLLABORATION

We are founded and led by committed members of the community of people living with disabilities and their caregivers. Five years and 26 patents later, our team continues to make this beautiful, dangerous world more accessible for all people using power wheelchairs. Along the way, we have created a culture of intentional partnership, because technology is another area where inclusion matters.

LUCI and the MyLUCI app are compatible with:

