

QUICKIE REQUEST FOR QUOTE



Please complete entire form

PROVIDER INFORMATION

Provider Company Name:	Billing Address:	Shipping Address:	Same as billing	
	RECIPIENT NAME:	RECIPIENT NAME:		
Provider Contact Name:	STREET:	STREET:		
Referring Therapist, Clinic and Contact Info:	CITY:	STATE:	CITY:	STATE:
Client Reference / Order Reference Number:	ZIP CODE:	COUNTRY:	ZIP CODE:	COUNTRY:

PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or a retrofit to an existing wheelchair?

Insurance (K0108 – Miscellaneous) or Private Pay?

New Chair Order

Retrofit Existing Chair

Insurance Carrier

Private Pay

Select LUCI Model

MSRP: \$9,854.00 USD

Note: LUCI is NOT compatible with the following Quickie power order form options:

- Captain Seat
- VR-2 Joystick
- No Power Seat Functions

Quickie Q700 M
LUCI Model 108
*We do not support 8 mph motors, Group 4 base, or seat depths <20".
Sedeo Pro req min 18.5" STFH.
All others req min 18" STFH.*

Quickie Q500 M
LUCI Model 107
Select Seat-To-Floor Height
Note: Measured from front of seat pan to floor with 0-degree angle.

Quickie Q300 M mini
LUCI Model 109

Select Special Modification:

Q700 UP Stander
Heavy Duty (HD) Model

HD Note: : LUCI is not designed to protect any part of the seat pan beyond the chair base.

16.5"	18"
17"	18.5"
17.5"	19"

Select Special Modification:

Zippie mini

Is this chair order from Sunrise Built4Me option?

Yes; Order #

If yes, please include the full chair order form and

No

LUCI team will confirm compatibility.

If you are using a configuration not found on this form, please call 615-813-5824, ext 1 before ordering.

Send completed form to orders@luci.com

SUNM-1224-C

QUICKIE REQUEST FOR QUOTE

Please complete entire form



Select Legrest + Footplate:

Note: If you intend to add a custom footbox, it should not exceed the width of the footplate. This will ensure adequate visibility for LUCI's camera sensors.

ADULT OPTIONS

Power Centermount, Standard Footplate (14" W x 10" D)

Power Centermount, Large Footplate (16" W x 12" D)

Fixed Centermount

Angle adjustable Centermount

Other:

Note: LUCI will confirm compatibility

PEDIATRIC OPTIONS

Fixed Centermount Legrest, Short (12.5" to 15.5" leg length)

Fixed Centermount Legrest, Long (14.5" to 19.5" leg length)

Power Centermount, Standard Footplate (14" W x 10" D)

Power Centermount, Large Footplate (16" W x 12" D)

Power Centermount Legrest (90° to 0°, 7.5" extension)

Angle Adjustable Centermount

Select Leg Length:

Small (3" to 7" leg length)

Medium (3" to 10" leg length)

Select Footplate:

Dual Small (5.75" x 6")

Dual Large (5.75" x 10")

Select Drive Method:

Note: This will ensure you receive the correct brackets

Standard LED Joystick

Color Joystick

Heavy Duty Color Joystick

High Resolution Color Joystick, (CJSM2+BT)

Omni2: ALT Drive Method

Note: We do not support Latch drive or VR2 devices

Select Power Seating Modules:

Note: We do not need to know if you have powered recline or legrests.

Seat elevator

Tilt

Tilt and seat elevator

No tilt or seat elevator

Any custom chair modifications? If so, explain.

REQUEST DELIVERY SUPPORT

LUCI is a finely-tuned, highly customizable system and works best when it is properly setup for the power wheelchair user.

Let LUCI be your guide during delivery and setup and schedule a 30-minute remote support session with us!

Scan here to schedule your session:

Delivery Contact Name:

Wheelchair User Name:



<https://calendly.com/d/2v7-77p-s64/luci-delivery-support>

If you are using a configuration not found on this form, please call 615-813-5824, ext 1 before ordering.

Send completed form to orders@luci.com

SUNM-1224-C