

QUICKIE REQUEST FOR QUOTE



Please complete entire form



CANADIAN PROVIDER INFORMATION

Provider Company Name:	Billing Address:	Shipping Address:	<i>Same as billing</i>	
	RECIPIENT NAME:	RECIPIENT NAME:		
Provider Contact Name:	STREET:	STREET:		
Referring Therapist, Clinic and Contact Info:	CITY:	PROV:	CITY:	PROV:
Client Reference / Order Reference Number:	POSTAL CODE:	COUNTRY:	POSTAL CODE:	COUNTRY:

PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or a retrofit to an existing wheelchair?

Insurance or Private Pay?

New Chair Order

Retrofit Existing Chair

Insurance Carrier

Private Pay

Select LUCI Model

MSRP: \$13,995.00 CAD

Note: LUCI is NOT compatible with the following Quickie power order form options:

- Captain Seat
- VR-2 Joystick
- No Power Seat Functions

Quickie Q700 M
LUCI Model 108

We do not support 8 mph motors, Group 4 base, or seat depths <20".
Sedeo Pro req min 18.5" STFH.
All others req min 18" STFH.

Select Special Modification:

Q700 UP Stander

Heavy Duty (HD) Model

HD Note: : LUCI is not designed to protect any part of the seat pan beyond the chair base.

Quickie Q500 M
LUCI Model 107

Select Seat-To-Floor Height
Note: Measured from front of seat pan to floor with 0-degree angle.

16.5"	18"
17"	18.5"
17.5"	19"

Quickie Q300 M mini
LUCI Model 109

Select Special Modification:

Zippie mini

Is this chair order from Sunrise Built4Me option?

If yes, please include the full chair order form and LUCI team will confirm compatibility.

Yes; Order #

No

If you are using a configuration not found on this form, please call 615-813-5824, ext 1 before ordering.

Send completed form to orders@luci.com

SUNM-0125-CAN

QUICKIE REQUEST FOR QUOTE



Please complete entire form



Select Footplate:

Note: If you intend to add a custom footbox, it should not exceed the width of the footplate. This will ensure adequate visibility for LUCI's camera sensors.

ADULT OPTIONS

- Power Centermount, Standard, 10" D x 14" W
- Power Centermount, Large, 12" D x 16" W
- Fixed Centermount
- Angle adjustable Centermount
- Other:

Note: LUCI will confirm compatibility

PEDIATRIC OPTIONS

- Fixed Centermount Legrest, Short (12.5" to 15.5" leg length)
- Fixed Centermount Legrest, Long (14.5" to 19.5" leg length)
- Power Centermount, Standard Footplate (14" W x 10" D)
- Power Centermount, Large Footplate (16" W x 12" D)
- Power Centermount Legrest (90° to 0°, 7.5" extension)
- Angle Adjustable Centermount
Select Leg Length: Select Footplate:
 - Small (3" to 7" leg length) Dual Small (5.75" x 6")
 - Medium (3" to 10" leg length) Dual Large (5.75" x 6")

REQUEST DELIVERY SUPPORT

LUCI is a finely-tuned, highly customizable system and works best when it is properly setup for the power wheelchair user.

Let LUCI be your guide during delivery and setup and schedule a 30-minute remote support session with us!

Select Drive Method:

Note: This will ensure you receive the correct brackets

- Standard LED Joystick
- Color Joystick
- Heavy Duty Color Joystick
- High Resolution Color Joystick, (CJSM2+BT)
- Omni2: ALT Drive Method

Note: We do not support Latch drive or VR2 devices

Select Power Seating Modules:

Note: We do not need to know if you have powered recline or legrests.

- Seat elevator
- Tilt
- Tilt and seat elevator
- No tilt or seat elevator

Any custom chair modifications? If so, explain.

Scan here to schedule your session:

Delivery Contact Name:

Wheelchair User Name:



<https://calendly.com/d/2v7-77p-s64/luci-delivery-support>

If you are using a configuration not found on this form, please call 615-813-5824, ext 1 before ordering.

Send completed form to orders@luci.com

SUNM-0125-CAN