

# AMYLOR REQUEST FOR QUOTE

Please complete entire form



## PROVIDER INFORMATION

Provider Company Name: Billing Address: Shipping Address: *Same as billing*  
RECIPIENT NAME: RECIPIENT NAME:  
  
Provider Contact Name: STREET: STREET:  
  
Referring Therapist, Clinic and Contact Info: CITY: STATE: CITY: STATE:  
  
Client Reference / Order Reference Number: ZIP CODE: COUNTRY: ZIP CODE: COUNTRY:

## PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or a retrofit to an existing wheelchair?

☐ New Chair Order

☐ Retrofit Existing Chair

If retrofitting a wheelchair, the serial number is required:

Select LUCI Model: [MSRP: \\$9,854.00 USD](#)

Note: We do not support HD models, 8.5 mph motors, 12-1/2" drive wheels, 9" front or rear casters, straight front caster arms, swing-away footrests, individual power elevating legrests, and extended battery compartments. Model requires at least one power seating function to have enough space to accommodate the LUCI hardware.

Amylor Alltrack R<sup>3</sup> Hybrid – LUCI Model 110

Select Foot Platform:

This will ensure you receive the correct brackets.

☐ Flip-up one piece Foot Platform

☐ Angle Adjustable for center mount/Flip-up Split Footplates for center mount

## REQUEST DELIVERY SUPPORT

LUCI is a finely-tuned, highly customizable system and works best when it is properly setup for the power wheelchair user.

Let LUCI be your guide during delivery and setup and schedule a 30-minute remote support session with us!

Select Drive Method:

This will ensure you receive the correct brackets

☐ R-net - (LED) – Joystick

☐ R-net Joystick with Color screen

☐ R-net - (LED) – Joystick with Lighting Operation

☐ R-net Joystick with Color screen and Lighting Operation

☐ R-net Advanced Bluetooth and Infrared Joystick (3.5" Screen)

☐ Omni2: ALT Drive Method

Note: We do not support Latch drive or VR2 devices

Any custom chair modifications? If so, explain.

Insurance (K0108 – Miscellaneous) or Private Pay?

☐ Insurance Carrier

☐ Private Pay

Scan here to schedule your session:

Delivery Contact Name:

Wheelchair User Name:



<https://calendly.com/d/2v7-77p-s64/luci-delivery-support>

Send completed form to [orders@luci.com](mailto:orders@luci.com)

If you are using a configuration not found on this form, please call 615-813-5824, ext 1 before ordering.

AMYL-0725-C