

QUICKIE REQUEST FOR QUOTE



Please complete entire form and send to orders@luci.com



CANADIAN PROVIDER INFORMATION

Provider Company Name:	Billing Address:	Shipping Address:	<i>Same as billing</i>	
	RECIPIENT NAME:	RECIPIENT NAME:		
Provider Contact Name:	STREET:	STREET:		
Referring Therapist, Clinic and Contact Info:	CITY:	PROV:	CITY:	PROV:
Client Reference / Order Reference Number:	POSTAL CODE:	COUNTRY:	POSTAL CODE:	COUNTRY:

PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or a retrofit to an existing wheelchair?

New Chair Order

Retrofit Existing Chair

Insurance or Private Pay?

Pro Tip: Visit luci.com/funding/ for insurance guidance resources, or connect with our team directly at education@luci.com.

Insurance Carrier

Private Pay

Select LUCI Model

MSRP: \$15,051 CAD

LUCI is ONLY compatible with R-Net Electronics. LUCI is NOT compatible with the following:

- Captain Seat
- Swing-Away Legrests
- VR-2 Joystick
- LiNX Electronics
- 8 mph motors
- Group 4 bases

Quickie Q700 M
LUCI Model 108
Max seat depth is 21"
Min STFH is 18.5"

Select Special Modification:

Q700 UP Stander

Heavy Duty (HD) Model*

**HD Note: LUCI is designed to protect what falls within the width of the chair's base. If the seat pan/arm rests fall outside of that width, they will not be protected by LUCI.*

Quickie Q500 M
LUCI Model 107

Select Seat-To-Floor Height
Measured from front of seat pan to floor with 0-degree tilt angle.

16.5"	18"
17"	18.5"
17.5"	19"

Quickie Q300 M mini
LUCI Model 109

Select Special Modification:

Zippie mini

Is this chair order from Sunrise Built4Me?

Yes; Order #

If yes, please include the full chair order form and LUCI team will confirm compatibility.

No

If you are using a configuration not found on this form, please call 615-813-LUCI before ordering.

Send completed form to orders@luci.com.

SUNM-0426-CAN

QUICKIE REQUEST FOR QUOTE

Please complete entire form and send to orders@luci.com



Select Footplate:

LUCI does not support swing-away footrests. This will ensure adequate visibility for LUCI's camera sensors.

ADULT OPTIONS

Power Centermount, Standard, 10" D x 14" W

Power Centermount, Large, 12" D x 16" W

Fixed Centermount

Angle adjustable Centermount

Other:

Note: LUCI will confirm compatibility

PEDIATRIC OPTIONS

Fixed Centermount Legrest, Short (12.5" to 15.5" leg length)

Fixed Centermount Legrest, Long (14.5" to 19.5" leg length)

Power Centermount, Standard Footplate (14" W x 10" D)

Power Centermount, Large Footplate (16" W x 12" D)

Power Centermount Legrest (90° to 0°, 7.5" extension)

Angle Adjustable Centermount

Select Leg Length:

Select Footplate:

Small (3" to 7" leg length)

Dual Small (5.75" x 6")

Medium (3" to 10" leg length)

Dual Large (5.75" x 6")

Select Drive Method:

This will ensure you receive the correct dashboard mounting bracket.

Standard LED Joystick

Color Joystick

Heavy Duty Color Joystick

High Resolution Color Joystick, (CJSM2+BT)

Omni2: ALT Drive Method

Note: We do not support Latch drive, VR2 devices, or LiNX Electronics.

Select Power Seating Modules:

We do not need to know if you have powered recline or legrests. Model requires at least one of the below power seating functions to have enough space to accommodate the LUCI hardware.

Seat elevator

Tilt

Will you be using lower extremity or foot support in addition to the standard footplates (e.g., foot box)?

Yes, please share details in the custom chair modifications box

No

Any custom chair modifications? If so, explain.

If you are using a configuration not found on this form, please call 615-813-LUCI before ordering.

Send completed form to orders@luci.com

SUNM-0426-CAN