



QUICKIE REQUEST FOR QUOTE

Please complete entire form and send to orders@luci.com



PROVIDER INFORMATION



U.S. Department of Veterans Affairs

Billing Address: Shipping Address: *Same as billing*
 RECIPIENT NAME: RECIPIENT NAME:

Ordering Clinic Name: STREET: STREET:

Ordering Clinician and Contact Info: CITY: STATE: CITY: STATE:

Client Reference / Order Reference Number: ZIP CODE: COUNTRY: ZIP CODE: COUNTRY:

PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or a retrofit to an existing wheelchair? New Chair Order Retrofit Existing Chair

<p>Select LUCI Model MSRP: \$9,845.00 USD</p> <p><i>LUCI is ONLY compatible with R-Net Electronics. LUCI is NOT compatible with the following:</i></p> <ul style="list-style-type: none"> • Captain Seat • Swing-Away Legrests • VR-2 Joystick • LiNX Electronics • 8 mph motors • Group 4 bases 	<p>Quickie Q700 M LUCI Model 108 <i>Max seat depth is 21"</i> <i>Min STFH is 18.5"</i></p> <p>Select Special Modification: Q700 UP Stander</p> <p>Heavy Duty (HD) Model*</p> <p><i>*HD Note: LUCI is designed to protect what falls within the width of the chair's base. If the seat pan/arm rests fall outside of that width, they will not be protected by LUCI.</i></p>	<p>Quickie Q500 M LUCI Model 107</p> <p>Select Seat-To-Floor Height <i>Measured from front of seat pan to floor with 0-degree tilt angle.</i></p> <table style="margin-left: auto; margin-right: auto;"> <tr> <td>16.5"</td> <td>18"</td> </tr> <tr> <td>17"</td> <td>18.5"</td> </tr> <tr> <td>17.5"</td> <td>19"</td> </tr> </table>	16.5"	18"	17"	18.5"	17.5"	19"	<p>Quickie Q300 M mini LUCI Model 109</p> <p>Select Special Modification: Zippie mini</p>
16.5"	18"								
17"	18.5"								
17.5"	19"								

Is this chair order from Sunrise Built4Me? Yes; Order #
 If yes, please include the full chair order form and LUCI team will confirm compatibility. No

If you are using a configuration not found on this form, please call 615-813-LUCI before ordering.

Send completed form to orders@luci.com.

QUICKIE REQUEST FOR QUOTE



Please complete entire form and send to orders@luci.com

Select Footplate:

LUCI does not support swing-away footrests. This will ensure adequate visibility for LUCI's camera sensors.

ADULT OPTIONS

Power Centermount, Standard, 10" D x 14" W

Power Centermount, Large, 12" D x 16" W

Fixed Centermount

Angle adjustable Centermount

Other:

Note: LUCI will confirm compatibility

PEDIATRIC OPTIONS

Fixed Centermount Legrest, Short (12.5" to 15.5" leg length)

Fixed Centermount Legrest, Long (14.5" to 19.5" leg length)

Power Centermount, Standard Footplate (14" W x 10" D)

Power Centermount, Large Footplate (16" W x 12" D)

Power Centermount Legrest (90° to 0°, 7.5" extension)

Angle Adjustable Centermount

Select Leg Length:

Select Footplate:

Small (3" to 7" leg length)

Dual Small (5.75" x 6")

Medium (3" to 10" leg length)

Dual Large (5.75" x 6")

Will you be using lower extremity or foot support in addition to the standard footplates (e.g., foot box)?

Yes, please share details in the custom chair modifications box

No

Any custom chair modifications? If so, explain.

Select Drive Method:

This will ensure you receive the correct dashboard mounting bracket.

Standard LED Joystick

Color Joystick

Heavy Duty Color Joystick

High Resolution Color Joystick, (CJSM2+BT)

Omni2: ALT Drive Method

Note: We do not support Latch drive, VR2 devices, or LiNX Electronics.

Select Power Seating Modules:

We do not need to know if you have powered recline or legrests. Model requires at least one of the below power seating functions to have enough space to accommodate the LUCI hardware.

Seat elevator

Tilt

If you are using a configuration not found on this form, please call 615-813-LUCI before ordering.

Send completed form to orders@luci.com.

SUNM-0426-V